

COST EFFECTIVE SYSTEM COVERING ALL STANDARD URODYNAMIC EXAMINATIONS



UNIQUE WIRELESS CONNECTION
AT 915 MHZ ENSURES STABILITY
AND NO SIGNAL DROPOUT

COMPATIBILITY WITH DIFFERENT TYPES
OF THE PRESSURE MEASURING SYSTEMS

Notify body No. 2265, 3EC International a.s.

MMTH_2

MEDKONSULT medical technology s.r.o.
Pasturova 15, 779 00 Olomouc, Czech Republic


T +420 581 113 030
E mmt@mmtsystems.com



Harmony

can work either as a small standalone uroynamics system or implemented as a wireless extension into the larger urodynamic system such as the Melody or Symphony.

The unit is equipped with 5 measuring channels with inputs for recording pressures, EMG and optional infusion weight.



The docking station for Harmony can be attached to either the Symphony cart or the Melody configuration.



GET INTO THE FLOW WITH MMT URODYNAMICS



Harmony

www.mmtsystems.com

www.mmtsystems.com

USER FRIENDLY SOFTWARE UDMvision

EASY TO INSTALL

PREDEFINED PROTOCOLS


METHODS CONFORM TO ICS STANDARDS

MULTIPLE NOMOGRAMS

CONFIGURABLE EXAMINATION METHODS
AND TEST REPORTS

PVR VALUE AUTOMATICALLY SAVED INTO TEST REPORTS

CONFIGURABLE LANGUAGES



WIRELESS CONNECTION WITH HOSPITAL NETWORK,
INCLUDING DICOM AND HL7

VOICE OR REMOTE CONTROL

ONLINE REMOTE SERVICE SUPPORT



Harmony
is the smallest mobile unit
for basic urodynamic examinations
with wireless connectivity.

MMT URODYNAMIC
SYSTEMS ARE
COMPATIBLE WITH
A WIDE RANGE
OF ACCESSORIES


Harmony

CYSTOMETRY



MMT infusion weight ensures
precise volume filling.

EMG / BIOFEEDBACK



EMG curve is smooth
and precise thanks
to an external preamplifier.

VALSALVA LEAK POINT PRESSURE*



Flexible camera allows easy
visualisation of urethral leaks.

*Third party product

3D BLADDER SCANNER*



The portable, intuitive, and accurate
bladder scanner provides real-time
PVR information. The PVR value
is available in seconds and saved
into test reports automatically.

*Third party product

UROFLOWMETRY



The waterproof flow meter
allows easy cleaning
and disinfection. The lithium
battery provides incredible
battery performance for
a truly wireless experience.
No need for daily charging.



Different types
of stand for flow meters
available - height
adjustable, designed to fit
underneath the commode
chair. Easy to assemble
folding version allows
easy travelling and more
flexibility.



Lightweight and easy to clean
commode chair for seated
voiding, suitable also
for videourodynamics
procedures.