




TECHNICAL SPECIFICATIONS

Flow sensor	Range	± 1.5 l/s
Metal mesh	Linearity	$< \pm 2\%$ in the range ± 1.5 l/s
Pressure	Range	± 2.5 kPa
	Linearity	1.5 %fs
Physical Dimensions	Height: 31 cm (12.2 in) Width: 29 cm (11.4 in) Depth: 26 cm (10.2 in)	
Weight	6.4 kg (14.7 lbs) with support arm. Device only 4.3 Kg (9.4 lbs)	
Testing signal mode	Single frequency: 5, 6, 8, 10 Hz Enhanced and optimized multiple frequency: 5 + 11 + 19 Hz Pseudo random noise (PSRN): 5-37 Hz	
Measurement accuracy	Rrs, Xrs: ± 0.1 cm H ₂ O/(l/s) or $\pm 10\%$ of the measured value Breathing pattern parameters: $\pm 10\%$ of the measured value Lung volumes: ± 50 ml or $\pm 2.5\%$ of the measured value Environmental conditions: Temperature: $\pm 1^\circ\text{C}$ Humidity: $\pm 3\%$ (Relative Humidity) Pressure: ± 100 Pa	
Calibration	Factory calibration according to international recommendations ERS (Eur Respir J. 2020;55(2):1900753). Daily calibration check with a test object (supplied with the device) and 3 liters calibration syringe (optional not supplied with the device), required for the measurement of slow spirometry volumes, closing volumes and airways specific conductance.	
Patient Load	0.68-0.95 cm H ₂ O/(l/s) in the frequencies of normal breathing (0.2-0.98 Hz)	
Dead Space	35 ml	
Display	10.1" HD color display with multi-touch capacitive touchscreen and anti-glare film (can be used with medical gloves)	
Electrical Specifications	Medical grade 100/240 V, 50/60 Hz 60 W input AC/15 VDC output power supply (included)	

Model: RESMON PRO FULL (REF: RT1100)

RESMON PRO FULL is distributed in your country by:

Some of the described features may not be approved for use in your country/region. Please contact your local distributor for further information. Measurement of SVC, IC, Closing Volume (CVtot) and specific airways conductance (sGrs) are not yet cleared by FDA for distribution in the USA.

Publications: The RESMON PRO FULL has been used in a wide selection of publications for clinical and research uses in asthma, COPD and evaluation of airways clearance therapies. Updated list available upon request.

RESTECH
RESPIRATORY TECHNOLOGY



Via M. Gioia 61/63
20124 Milano, Italy

P.IVA/CF: 07116980967
REA: MI-1936811



www.restech.it

RESMON PRO FULL



ENCOURAGE ADOPTION
FOR BETTER RESPIRATORY DIAGNOSIS
MAKE
OSCILLOMETRY
ROUTINE
RESTECH

RESTECH
As simple as breathing



INDICATION FOR USE

THE RESMON PRO FULL is intended to measure respiratory system impedance using the Forced Oscillation Technique (FOT). RESMON PRO FULL is intended for use with pediatric and adult patients 3 years of age or older. The device is designed to be used by pulmonologists, general practitioners, nurses, respiratory therapists, laboratory technologists, medical researchers and similarly trained personnel in hospitals, clinics, and private physician offices.

ACCURATE AND NON-INVASIVE


RESMON PRO FULL utilizes oscillometry (also known as forced oscillation technique or FOT) to accurately and non-invasively measure the mechanical properties of the respiratory system using just a few normal tidal breaths.

This advanced technology offers a simple and effort-independent method to detect and localize airway obstruction, allowing differentiation between peripheral, central, or heterogeneous types.

ADVANCED FEATURES FOR ENHANCED PERFORMANCE

RESMON PRO FULL is equipped with advanced calculation algorithms based on the latest scientific research and clinical findings, ensuring exceptional performance.

By integrating respiratory physiology with state-of-the-art technology, RESMON PRO FULL introduces a range of new advanced parameters that combine oscillometry and static lung volume measurements, such as closing volume (CV_{fo}), critical reactance (X_{crit}), and oscillometric airway conductance (sGr_s). This comprehensive set of measurements enhances the physician's ability to detect diseases, even in their earliest stages, and effectively monitor the impact of therapeutic interventions on disease progression. Measurement of Closing Volume (CV_{fo}), critical reactance (X_{crit}) and oscillometric airway conductance (sGr_s) are not yet cleared by FDA for distribution in the USA.



DYNAMIC “WITHIN-BREATH” TESTING

- Real-time display of resistance (Rrs), reactance (Xrs) and flow or volume (user selectable)
- Measurement of inspiratory, expiratory and total Rrs and Xrs parameters
- Spectral analysis with the calculation of R5-19, AX and Fres
- Automatic detection and quantification of tidal expiratory flow limitation (EFLt) with ΔXrs index graph and % of flow limited breaths (FL%), (patent nr. wo2003103493)
- Full respiratory pattern measurement and reporting of VE, Vt, RR, Ti/Ttot, Vt/Ti, Vt/Te
- Customizable tidal loops graphs (available signals: Flow, Volume, Resistance, Reactance and Impedance)

SLOW VITAL CAPACITY TESTING

- SVC for monitoring of restrictive patterns
- IC for hyperinflation detection, pre-post testing for evaluating effect treatment
- Closing volume (CV_{fo}) and critical reactance (X_{crit}) for the evaluation of de-recruitment of peripheral airways
- Specific airway conductance (sGr_s) for the evaluation of therapeutic intervention compensated for changes in lung volumes
- Measurement of SVC, IC, Closing Volume (CV_{fo}) and specific airways conductance (sGr_s) are not yet cleared by FDA for distribution in the USA.

THREE MEASUREMENT MODES

- Enhanced and optimized multi-frequency mode of 5-11-19 Hz, from children to adults
- Single-frequency mode options of 5, 6, 8, 10 Hz (for children, severely obstructed patients and special research purposes)
- Enhanced and optimized Pseudo Random Noise (PSN) of 5-37 Hz

ENHANCED COMFORT AND TESTING EFFICIENCY

- Automatic adjustment of stimulus amplitude based on patient impedance to maximize measurement quality and optimal patient comfort
- Efficient testing with automatic discard of breaths affected by artifacts (i.e. cough, glottis closure, etc) and choice of number of accepted breaths, or time to end the test automatically, minimizing operator intervention
- CoV % (coefficient of variability) within-measurement and within-session with color coded warning, for optimal quality control
- Operator accessible screen, post-test, to review and edit breaths to include into the calculation of the results
- User selectable sets of reference equations for predicted normal value and Z-score calculation

CONVENIENT SIZE AND FUNCTIONS

- Self-contained and compact unit with data transfer capability on USB memory stick, directly to PCs via USB-OTG or LAN cable
- Connectivity options, USB and HDMI ports, with external keyboard, mouse, screen or projector
- Active-wash fan for expired CO₂ removal from the system reducing system dead space to 35 mL
- Unique impedance verification tube included, to daily check the accuracy of both resistance (Rrs) and reactance (Xrs)
- Wide 10.1" color touchscreen for fast, easy and intuitive test management
- Built-in database to store, retrieve, edit or delete patient data
- Multi-user access with password protected data security
- Detailed report of results for trending, analysis and statistics. Data available in raw format, XSL/CSV and PDF output as well as export to external PC (via USB-OTG or LAN cable for safe transfer connection and data protection)

PRE-POST RESULTS

