



Minispir[®]

Computer-based
Spirometer



Real time Flow/Volume and Volume/Time curves with PRE/POST Bronchodilator comparison on your PC using Winspiro PRO[®] software

Includes a wide range of selectable parameters

Spirometry test interpretation

Embedded Temperature sensor for BTPS conversion



Always included: Winspiro PRO[®] PC software with free update



Available with both DISPOSABLE or REUSABLE turbine flowmeter



Spirometry test: FVC, VC, IVC, MVV, PRE/POST Bronchodilator comparison with a wide range of selectable parameters

MINISPIR[®] PARAMETERS

Spirometry parameters

FVC, FEV1, FEV1/FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, Lung Age, Extrap. Volume, FET, FEV3, FEV3/FVC, FEV6, FEV1/FEV6, FIVC, FIV1, FIV1/FIVC, PIF, VC, IVC, IC, ERV, Rf, VE, VT, t_i, t_E, VT/t_i, t_E/t_{TO}, MVV

Minispir Light[®]

Simplified Computer-based
Spirometer



Real time Flow/Volume and Volume/Time curves on your PC using Winspiro Light[®] software

COPD and Asthma screening have never been so intuitive and inexpensive

Includes a wide range of selectable parameters

Spirometry test interpretation

Embedded Temperature sensor for BTPS conversion



Always included: Winspiro Light[®] simplified PC software with free update



Function available: PRE/POST Bronchodilator software feature



Available with DISPOSABLE turbine flowmeter



Spirometry test: FVC, VC, IVC

MINISPIR LIGHT[®] PARAMETERS

Spirometry parameters

FVC, FEV1, FEV1%, FEV6, PEF, FEF25-75, FIVC, Lung Age, VC, IVC



Ph: 1300 009 663 Fax: 1300 099 300

Email: info@zonemedical.com.au

7/22 Mavis Court, Yatala, Qld, 4207

www.zonemedical.com.au