

BedSide Infotainment Terminal

FUNCTION

- Plug-in module with multiple parameters available.
- Multiple display modes including other beds view and big font.
- Optional 12 or 15-inch screen.
- WiFi connectivity for flexible
- application.
- eTeleViewTM Central Monitor
- · system support.
- Smart alarm system for enhanced safety and efficiency.
- Fanless design for quiet operation.

DESCRIPTION

Multi-scenario Application:

Animal Care

- · Operating room use
- Plug-in Modular system

Central Monitoring System:

- Supports up to 64 beds or patients simultaneously.
- 72 hours of 64-channel holographic physiological waveform storage and review.
- Stores up to 240 hours of trend data and 720 alarm events per bed.
- Bi-directional communication with W12M/W15M Series monitors for enhanced patient care.







- Shortcut Menu
- Notification Panel
- Application Interface

Power On/Off:

- Power On
- Power Off

DEVICE FUNCTION

Usage Of Settings:

- Application Settings
- Display Settings
- Sound Settings
- · Security & Location
- Accessibility
- Digital Wellbeing & Parental Controls
- System
- Using COM Port
- Using NFC/IC
- Root Access Switch Steps
- · Device Files Transfer
- SD Card Usage Instruction

Standard Accessories

- DC 12V power adaptor
- · User manual card



