




BARIATRIC Evacuation

Bariatric and ICU patients can pose an especially difficult problem in the event of an evacuation. Help your staff by having the proper equipment available. 10% of your evacuation devices should be bariatric capable to meet minimum requirements.

#MS48
Bariatric Med Sled®
 Deployed Dimensions:
 48"W x 87"L
 Stored Dimensions:
 9" dia. x 48"W



#MS36
Standard Med Sled.
 Deployed Dimensions: 36"W x 87"L
 Stored Dimensions: 9" dia. x 36"W

#MS48
Bariatric Med Sled.
 Deployed Dimensions: 48"W x 87"L
 Stored Dimensions: 9" dia. x 48"W

KEY BENEFITS

- + Holds 1,000 lbs.
- + 48" wide to accommodate a patient's wide girth
- + Effective for ICU and Critical Care to add extra cushion and critical equipment
- + Typically used for patients of 300+ lbs.



#MSSDB2 **2-Sled Bariatric Storage Device**
 16"W x 48"H x 8 1/2"D

- + **Accommodates:** (2) Bariatric MS48 Sleds® positioned side-by-side (alternate option: 1 Standard MS36, 1 Bariatric MS48)

#MSSDB6 **6-Sled Bariatric Storage Device (Not shown)**
 50"W x 48"H x 8 1/2"D

- + **Accommodates:** (6) Bariatric MS48 Sleds® positioned side-by-side

EASY TO STORE
 WALL UNIT STORAGE DEVICES
 AVAILABLE FOR SLEDS
 (Sleds sold separately, not included
 with storage bag)

ORDERING INFORMATION

+ **Sales and Customer Service:** 1-866-207-5993 or for more information:
 Sales@medsled.com