



Value-Lift Treatment

A versatile table suitable for a wide range of modalities, including light physiotherapy use, massage, acupuncture, osteopathy, beauty treatments and examinations. The three-section configuration features a 700mm backrest and a long seat section. Ideal for treatments that require extra leg room.



GAS-ASSISTED ADJUSTABLE BACKREST

Lift-assist gas springs help support a patient's body weight for more comfortable and secure table adjustments. Patients can enjoy smoother table transitions with minimal disruption to treatment flow.



GAS-ASSISTED ADJUSTABLE HEADREST

The standard configuration includes a face hole in the head section for comfortable prone treatments. The head section can be finely adjusted to enable patients to relax neck muscles fully during treatments.



STAINLESS STEEL FOOTSWITCH

The table comes with a sturdy footswitch made with a stainless steel body. Order an additional footswitch to control the height from either side of the table.



DANISH DESIGNED LINAK ACTUATOR

All Value-Lift tables feature a Danish actuator by Linak which is manufactured to our strict quality control. All products are inspected and checked by Athlegen when they arrive in Australia.



ADJUSTABLE FLOOR-LEVELLING FEET

The table sits securely on four solid supports that can be adjusted for an uneven floor.



RETRACTABLE WHEELS

The foot operated retractable wheel system allows for easy manoeuvrability when moving the table. When not in use the wheels retract to provide a stable solid base.

KEY SPECIFICATIONS

CUSHIONING	60mm, medium finish
UPHOLSTERY	Semi-PU fabric
ADJUSTMENTS	Backrest: +65° Head section: +30° to -30°
SAFE WORKING LOAD	165kg
MOTOR	Origin: PRC Design: Linak Denmark Lift force: 6,000N Speed: 12mm/sec Height adjustment (min to max): 295mm in 26 seconds Speed test load: 165kg
FOOTSWITCH	Single or optional dual footswitches (one on each side of the table) Origin: PRC Stainless steel Power control
SHIPPING	1920 x 750 x 720mm 74kg
WARRANTY	Structural: 5 years Hardware: 3 years Motor: 2 years Control box: 2 years Footswitch: 2 years Gas strut: 2 years Upholstery: To arrive in perfect condition



