



## REHAB TILT IN SPACE MOBILE SHOWER COMMUNE WITH SWINGAWAY FOOTREST - 500MM

Item No. KA11512-REHAB

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind . Side access standard on 445 and 500mm seat widths . Lateral thoracic support mounting plates . 10-35° tilt angle . Hip/thigh guide support mounting plates . Height adjustable frame . Tension



## KCare Healthcare Solutions - Tilt In Space Rehab Shower Commode



### DESCRIPTION

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind .
- Side access standard on 445 and 500mm seat widths .
- Lateral thoracic support mounting plates .
- 10-35° tilt angle .
- Hip/thigh guide support mounting plates .
- Height adjustable frame .
- Tension adjustable backrest .
- Integrated push-handle for manoeuvrability .
- Height adjustable flip up arms for unrestricted side transfer .
- Swingaway and removable height adjustable footrest with contoured padded calf supports .
- 4x 125mm total locking castors .
- Designed to be used with a K Care patented removable seats .
- Available in 400, 445, 500mm seat widths .
- Padded seat not included .

## SPECIFICATIONS

Overall Width	635 mm
Overall Depth	1070 mm
Seat Width	500 mm
Seat Depth	435 mm
Min Floor To Seat Height	465 mm
Max Floor To Seat Height	540 mm
Armrest Type	Loop
Backrest Height	855 mm
Footrest Type	Swing Away
Product Weight	23 kg
Max User Weight	200 kg
Brand	K Care