

SYS-3010 Syringe Pump

Easy to use

Smart Pump with Clear and Simple key pad design, LCD Screen, friendly interface and complete software and hardware design. Users can operate it very easily even without the instruction of manual.

Reliable

SYS-3010 Syringe Pump is based on the IEC/ISO standard and in accordance with the rules of the FDA guide. It has a patent accurate transmission mechanism, more internal monitoring sensors, and a waterproof mechanism.





Easy operation and safety

- Friendly interface
- Stable drive mechanism with patent
- Accurate 16-step motor
- Up to 11-sensor smart monitoring
- Up to 13-point waterproof design
- Automatic pressure decrease during tube block







Dual independent sound-light alarm system

Dual CPU control system, better control and data transmission

Multi function interface, AC/DC input, wire/wireless data connection, nurse call

Main parameter

Power supply	AC: 100-240V; external DC: 12V; built-in lithium battery, at least 8 hours
Compatible syringe	10ml, 20ml, 30ml and 50 (60) ml
Infusion speed range	0.1-2000mL/h, 0.1mL/h minimum increment
Infusion mode	Rate, Time, Bodyweight Mode
Infusion accuracy	Mechanism accuracy: ±1%, Accuracy with syringe: ±2%
Bolus rate	300-1500ml/h
KVO	0.1-2.5ml/h

MEDCAPTAIN[]

MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.



Add: 12th Floor, Baiwang Research Building, No.5158 Shahe West Road, Xili, Nanshan, 518055 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA

Tel: +86 755 26953369 E-mail: info@medcaptain.com Fax: +86 755 26001651 Website: www.medcaptain.com

MEDCAPTAIN is a trade mark of Medcaptain Medical Technology Co., Ltd. Specifications are subject to change without prior notice.

©2017-2018 Medcaptain Medical Technology Co., Ltd. All rights reserved.

ENG-SYS-3010-2P-20170427





