

Portable Oxygen Concentrators

zone
M E D I C A L


1300 009 663

info@zonemedical.com.au

If you, or someone you're caring for has been prescribed supplemental oxygen therapy, Zone Medical is here to assist you in finding the best solution tailored to your needs. We offer a variety of lightweight portable oxygen concentrators, perfect for individuals on the go or traveling. Our battery operated models provide the freedom of movement, allowing you to embrace life to the fullest.

Zone Medical's portable oxygen concentrators are available in continuous or pulse dose models. Pulse dose models deliver a gentle shot of oxygen upon inhalation and feature multiple settings to adjust the size of the dosage. It's essential to consult with your healthcare provider to determine the appropriate setting for your device.

When traveling by air, the altitude can impact your oxygen saturation levels, potentially necessitating a higher setting on your device. Discuss with your healthcare provider the appropriate setting for flying. We recommend contacting your airline in advance to obtain approval and secure a seat with a power outlet for using your device during the flight.



1000ml/minute on pulse delivery is equivalent to approx 5 LPM continuous flow delivery.

All Models are FAA approved for domestic and international flight.



Model Name	Kington P2	Inogen Rove 6	Philips Simply Go Mini	Caire Freestyle Comfort
Brief Description	Smallest and lightest in the world. Includes very helpful additional accessories as standard. Best bang for your buck.	Most powerful and longest battery life for a pulse unit under 3kg. The original and best seller globally.	Portable and sleek design, offering pulse dose oxygen delivery. Lightweight and powerful.	Ergonomically designed for comfort and ease of use. The curvature of the machine fits nicely on the body.
Delivery Mode	Pulse Dose	Pulse Dose	Pulse Dose	Pulse Dose
Pulse Dose Settings	5 Settings	6 Settings	5 Settings	5 Settings
Max Pulse Dose Per Minute (mL)	1000mL	1260mL	1000mL	1050mL
Continuous Flow Rate (LPM)	-	-	-	-
Size With Battery (HxWxD cm)	16.0 x 22.1 x D8.5	20.5 x 18.3 x 8.3	23.9 x 21.1 x 9.1	25.4 x 18.5 x 7.9
Weight With Battery	1.98kg	2.2kg (single battery) 2.6kg (double battery)	2.3kg (single battery) 2.7kg (double battery)	2.3kg (single battery) 2.7kg (double battery)
Noise Decibel	≤ 42dB	≤ 39dB	≤ 43dB	≤ 43dB
Battery Run Time (on Setting 2)	~ 4hrs	Up to 6.25hrs (single battery) Up to 12.75hrs (double battery)	Up to 4.5hrs (standard battery) Up to 9hrs (extended battery)	Up to 4 hours (single battery) Up to 8 hours (double battery)
Recharge Time	4 hrs via Machine 2 hrs via External Battery Charger	Up to 3hrs (single battery) Up to 6hrs (double battery)	Up to 4hrs (standard battery) Up to 8hrs (extended battery)	Up to 3.5hrs (single battery) Up to 6hrs (double battery)
Included Accessories	2 Batteries, Carry Bag, AC Charger, DC Charger, External Battery Charger, 10 Nasal Cannula, Pulse Oximeter, backpack	1 Battery, Carry Bag, AC Charger, DC Charger	1 Battery, Carry Bag, Accessory Bag, AC Charger, DC Charger	1 Battery, Carry Bag, AC Charger, DC Charger
Warranty (may change depending on distributor)	Device - 3 Years Battery - 1 Year Accessories - 1 Year Sieve Beds - 1 Year	Device - 3 Years Battery - 1 Year Accessories - 1 Year Sieve Beds - 1 Year	Device - 3 Years Batteries - 1 Year Accessories - 1 Year Sieve Beds - 1 Year	Device - 3 Years Battery - 1 Year Accessories - 1 Year Sieve Beds - 3 Years

Note: single battery and 8-cell, and double battery and 16-cell, are sometimes used interchangeably. Wording used is dependent on the manufacturer.



Model	GCE Zen-O Lite	Caire Eclipse 5	Philips SimplyGo	GCE Zen-O
Brief Description	Lightweight and portable, designed for maximum mobility. Provides pulse dose oxygen delivery, making it ideal for active users on the go.	Most powerful unit on the market for both continuous flow and pulse dose delivery. Ideal for home and travel, with a comprehensive range of settings.	Compact and durable with both continuous flow and pulse dose options. Known for its reliability and ease of transport, making it perfect for everyday use.	Dual-mode operation with both continuous and pulse flow settings. Sturdy build and excellent battery life, offering flexibility and independence for users
Delivery Mode	Pulse Dose	Pulse Dose / Continuous	Pulse Dose / Continuous	Pulse Dose / Continuous
Pulse Dose Settings	5 Settings	9 Settings	5 Settings	6 Settings
Max Pulse Dose Per Minute (mL)	1050mL	2880mL	1260mL	1980mL
Continuous Flow Rate (LPM)	-	0.5 - 3L / min	0.5 - 2L / min	0.5 - 2L / min
Size With Battery (HxWxD cm)	24.9 x 23.1 x 9.7	49.0 x 31.2 x 18.0	29.2 x 25.4 x 15.2	21.2 x 31.3 x 16.8
Weight With Battery	2.5kg	8.3kg	4.54kg	4.66kg
Noise Decibel	≤ 37dB	≤ 40dB	≤ 43dB	≤ 38dB
Battery Run Time (on Setting 2)	Up to 4hrs	Up to 5.1hrs (pulse dose setting 2) Up to 2hrs (continuous flow)	Up to 3hrs (pulse dose setting 2) Up to 54min (continuous flow)	Up to 4 - 8hrs* (pulse dose setting 2) Up to 0.75 - 1.5hrs* (continuous flow)
Recharge Time	Up to 3hrs	~ 2 - 5hrs (while operating)	~ 2 - 3hrs (while operating)	~ 3 - 6hrs* (while operating)
Included Accessories	1 Battery, Carry Bag, AC Charger, DC Charger	1 Battery, Trolley, AC Charger, DC Charger	1 Battery, Carry Bag, AC Charger, DC Charger, Accessory bag	2 Batteries, Carry Bag, Trolley AC Charger, DC Charger
Warranty (may change depending on distributor)	Device - 3 Years Battery - 1 Year Accessories - 1 Year Sieve Beds - 1 Year	Device - 3 Years Batteries - 1 Year Accessories - 1 Year Sieve Beds - 1 Year	Device - 2 Years Batteries - 3 Months Accessories - 3 Months Sieve Beds - 1 Year	Device - 3 Years or 15,000hrs Batteries - 1 Year Accessories - 1 Year Sieve Beds - 1 Year

* Depending on type of battery installed.

Variation depending on the type of battery installed. Single battery (or 8-cell) vs double battery (or 16-cell).