

## Product Recommendation at the Point of Sale







- Measurement of skin hydation oil content skin elasticity melanin content & need of sun protection products
- Display & interpretation of the results on diode chains
- Visioscope® skin camera for impressive close up images of skin and hair on the monitor

## Visioface 1000 D

- Very homogeneous illumination to take images of the face
- High resolution camera to see all details of the skin
- Comprehensive, economic analysis of the face: wrinkles, spots, pores, colour differences
- UV-like image for hidden spots



- Flexible customized solutions of hard- and software to suit the needs of your brand
- From small handheld battery powered devices to smart multiparametric systems with sophisticated software
- Our worldwide acknowledged skin testing instruments help you get ahead of your competitors





