



# REHAB TRANSIT MOBILE SHOWER COMMUNE WITH SWINGAWAY FOOTREST - 320MM

Item No. KA114SC-REHAB

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind . Side access standard on 445 and 500mm seat widths . Lateral thoracic support mounting plates . Hip/thigh guide support mounting plates . Height adjustable frame . Tension adjustable backrest .



## DESCRIPTION

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind .
- Side access standard on 445 and 500mm seat widths .
- Lateral thoracic support mounting plates .
- Hip/thigh guide support mounting plates .
- Height adjustable frame .
- Tension adjustable backrest .
- Integrated push-handle for manoeuvrability .
- Height adjustable flip up arms for unrestricted side transfer .
- Swingaway and removable height adjustable footrest .
- 4x 125mm total locking castors .

- Designed to be used with a K Care patented removable seats .
- Available in 320, 400, 445, 500mm seat widths .
- Padded seat not included .

## SPECIFICATIONS

Overall Width	485 mm
Overall Depth	750 mm
Seat Width	320 mm
Seat Depth	320 mm
Min Floor To Seat Height Front	475 mm
Max Floor To Seat Height Front	575 mm
Min Floor To Seat Height Rear	452 mm
Max Floor To Seat Height Rear	552 mm
Armrest Type	Loop
Backrest Height	415 mm
Footrest Type	Swing Away
Product Weight	12 kg
Max User Weight	120 kg
Brand	K Care